

<u>Pre Clean Flow Test</u>	<u>Post Clean Flow Test</u>	<u>Pre Clean Flow Test cc/min</u>	<u>Post Clean Flow Test cc/min</u>
133.5		534	0
118.5		474	0
119.25		477	0
124.5		498	0
141		564	0
133.5		534	0
		0	0
		0	0

Date: 12/30/2008
Customer Number: 4073835118
Sold To:

Attention:


	<u>Leak Test</u>	<u>Resistance</u>	<u>Spray Pattern</u>	<u>Pre Clean Flow Test</u>	<u>Post Clean Flow Test</u>	<u>Flow Increase</u>
Injector 1	Ok	13	Fair	534	0	-100.00%
Injector 2	Ok	12.6	Good	474	0	-100.00%
Injector 3	Ok	12.6	Good	477	0	-100.00%
Injector 4	Ok	15.5	Fair	498	0	-100.00%
Injector 5	Ok	13	Good	564	0	-100.00%
Injector 6	Ok	12.6	Good	534	0	-100.00%
Injector 7				0	0	-100.00%
Injector 8				0	0	-100.00%

Injector Notes:

1=JW, 2=S08, 3=S09, 4=RR, 5=RD,
6=Black
Left picture is JW, S08, S09...right picture is RR injector

Spray Pattern Pre-Cleaning

Spray Pattern Post-Cleaning



Injector parts replaced:

- Pintle Cap
 Filter
 Top O-Ring
 Bottom O-Ring
 Spacer
 Head Seal
 Clamp Seal
 Hose